

ELECTRIC HEATING INCUBATOR



BPX series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

ELECTRIC HEATING INCUBATOR >>

BPX series

Product Features >

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted optionally;
- Adopting PMMA I operating system, temperature control, timing, and over temperature alarm;
- Imported motor and fan, intelligent temperature controller with timing function;
- The fan has 6-stage speed regulation;
- LCD screen, temperature and time control, over-temperature alarm function;
- Observation window on the door, convenient for users to observe;
- Nano material magnetic door seal strips and heat preservation material;
- Programmable design, it can be set 10 segments 100 cycles;
- Power supply recovery function, data will not be lost when power failure and system halted.



Shelf



Temperature Sensor



Control Panel



Model	BPX-52	BPX-82	BPX-162	BPX-272	BPX-452	BPX-652	BPX-1002
Volume [L]	40	76	142	234	452	650	1000
Display	LCD						
Temp. Range [°C]	RT+5~85				RT+5~100		
Temp. Resolution [°C]	0.1						
Temp. Fluctuation at 37°C [°C]	±0.5				±0.8		
Temp. Uniformity at 37°C [°C]	±1				±2		
Power Supply [V, Hz]	220±10%, 50±2%						
Power [W]	300	450	450	750	1500	1500	2000
Timer Range (min/h)	0~9999						
Internal Dimensions W×D×H [mm]	350×330×350	450×380×450	550×470×550	600×520×750	600×580×1300	720×650×1400	1100×650×1400
Housing Dimensions W×D×H [mm]	640×556×505	740×606×605	840×696×705	890×746×905	774×845×1704	865×866×1730	1220×850×1730
Number of Shelves (std./max.) [pc]	2/4	2/5	3/7	3/10	4/18	4/18	4/18
Net Weight /Gross Weight [kg]	37/47	50/61	68/91	88/111	-	-	-
Package Size [mm]	740*696*645	840*746*745	960*856*855	1010*906*1055	-	-	-