

WATER-JACKET



BG Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

WATER-JACKET

BG Series

Product Features

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted optionally;
- LCD screen, temperature and time control, over-temperature alarm function;
- Heating with water jacket at both sides and bottom of the working chamber;
- A high quality low-noise fan is equipped in the working chamber to ensure the temperature uniformity;
- Power supply recovery function, data will not be lost when power failure and system halted.

Model	BG-50	BG-80	BG-160	BG-270
Volume [L]	46	70	160	247
Temp. Range [°C]	RT+2-65			
Temp. Resolution [°C]	0.1			
Temp. Fluctuation at 60°C [°C]	±0.2			
Power Supply [V, Hz]	220±10%, 50±2%			
Power [W]	430	630	1030	1530
Internal Dimensions W×D×H [mm]	310×380×450	400×400×500	550×500×650	600×600×750
Housing Dimensions W×D×H [mm]	460×510×765	550×530×815	700×630×965	750×730×1065
Number of Shelves (std./max.) [pc]	2/4	2/4	3/5	3/6
Net Weight /Gross Weight [kg]	40/48	49/58	80/98	95/115
Package Size [mm]	580*620*910	670*640*960	830*750*1120	890*860*1260