



BIOCHEMICAL INCUBATOR



BSP Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

BIOCHEMICAL INCUBATOR >>

BSP Series

Product Features

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted optionally;
- PMMA II operating system, with over temperature alarm and self-diagnosis function;
- Color LCD backlight display
- With screen lock time, screen saver password setting function, to prevent operation at will;
- Programmable design, it can be set 30 segments 99 cycles;
- Use centrifugal fan to make the temperature of the working room uniform and stable;
- Observation glass inner door;
- Test holes are equipped on both sides, internal diameter is 30mm;
- The fan has 6-stage speed regulation,
- With reservation and timing function;
- Equipped with mechanical lock;
- Equipped with 220V power socket;
- Equipped with fluorescent lamp in working chamber (mould incubator is equipped with UV light);
- Adopt imported compressor unit, low noise, stable performance;
- Power supply recovery function, data will not be lost when power failure and system halted.



Shelf



Inner Glass Door Lock



Glass Lock Group

Model	BSP-100	BSP-150	BSP-250	BSP-400	BSP-800
Volume [L]	100	150	250	400	800
Type	Biochemical Incubator				
Temp. Range [°C]	4~70				
Temp. Resolution [°C]	0.1				
Temp. Fluctuation [°C]	±0.5°C (heating) ±0.7°C (cooling)				
Temp. Uniformity at 37°C[°C]	±1			±1.5	
Power [W]	1000	1000	1100	1350	2000
Timer Range [min/h]	0~9999				
Reservation Range [min]	0~9999				
Internal Dimensions W×D×H [mm]	490×390×610	510×390×760	510×450×1090	600×640×1050	1220×585×1123
Housing Dimensions W×D×H [mm]	630×680×1250	650×680×1400	650×740×1726	745×930×1695	1503×910×1820
Number of Shelves (std./max.) [pc]	2 / 4	3 / 5	4 / 7	4 / 7	6 / 12
Load per Shelf [kg]	15				
Net Weight/Gross Weight [kg]	91/109	103/119	120/150	150/182	286/342
Package Size [mm]	730*820*1440	750*820*1590	750*880*1910	890*1120*1880	1640*1090*2040